

BGI Sample Testing Report

1. Project Information

Report No.: THKa15121103

Project Name	UW - Sam **** - Sam White - de-novo genome Ostrea lurida	Project No.	F15FTSUSAT0327			
Customer Name	Steven Roberts	Customer Unit	Univ. of Washington			
Lab Sample Collector	WONGWINGLEUNG	Lab Sample Receiving Date	20151211			
Lab Sample Tester	yutszfung, wongwingleung	Lab Sample Testing Date	20151211			
Reported by	yutszfung Inspected yutszfung	Approved by	Wong Wai Man Report Vivien Date 20151211			

2. Sample Test Method

- ①Method of concentration determination:■Qubit Fluorometer、□NanoDrop、□Microplate Reader
- ②Method of sample integrity test: ■Agarose Gel Electrophoresis

3. Sample Test Result

No.	Sample Name	Sample Number	Tube No.	Concen- tration(ng/ μL)	Volume(μL)	Total Mass(µg)	sample Integrity	Library type	Test result	Remark
1	oyster_01	8521512001643	1	154	160	24. 64	Degraded Moderate	10K Mate Pair	Level D	The sample is degraded moderately. the total mass is too low. Proposed to resend the sample.
2	oyster_01	8521512001643	1	154	160	24. 64	Degraded slightly	5k/6K Mate Pair	Level B	
3	oyster_01	8521512001643	1	154	160	24. 64	Degraded slightly	2K Mate Pair	Level B	
4	oyster_01	8521512001643	1	154	160	24. 64	Degraded slightly	PCR-free	Level A	The sample is degraded slightly.
5	oyster_01	8521512001643	1	154	160	24. 64	Degraded slightly	≤800bp Insert Size	Level A	
6	oyster_02	8521512001644	1	428	150	64. 2	Degraded slightly	≤800bp Insert Size	Level A	
7	oyster_02	8521512001644	1	428	150	64. 2	Degraded slightly	PCR-free	Level A	The sample is degraded slightly.
8	oyster_02	8521512001644	1	428	150	64. 2	Degraded slightly	2K Mate Pair	Level A	
9	oyster_02	8521512001644	1	428	150	64. 2	Degraded Moderate	10K Mate Pair	Level C	The sample is degraded moderately. Proposed to resend the sample.
10	oyster_02	8521512001644	1	428	150	64. 2	Degraded slightly	5k/6K Mate Pair	Level A	
11	oyster_03	8521512001645	1	296	53	15. 69	Degraded slightly	PCR-free	Level B	The sample is degraded slightly.
12	oyster_03	8521512001645	1	296	53	15. 69	Degraded slightly	2K Mate Pair	Level C	The total mass is too low. Proposed to resend the sample.
13	oyster_03	8521512001645	1	296	53	15. 69	Degraded slightly	5k/6K Mate Pair	Level D	The total mass is too low. Proposed to resend the sample.
										The sample is degraded

14	oyster_03	8521512001645	1	296	53	15. 69	Degraded Moderate	10K Mate Pair	Level D	moderately. the total mass is too low. Proposed to resend the sample.
15	oyster_03	8521512001645	1	296	53	15. 69	Degraded slightly	≤800bp Insert Size	Level A	

Note*:

- 1. The test result based on the 《DNA sequencing sample quality standards》 explains whether the testing sample meets the requirement of library construction.
- a) Level A means the sample is qualified, and the amount of sample is sufficient for two or more library constructions.
- b) Level B means the sample is qualified, but the amount of sample only satisfies one time library construction.
- c) Level C means the sample does not totally meet the requirements of library construction and sequencing. BGI can try to construct the library but the quality of the sequence is not guaranteed
- d) Level D means the sample does not meet the requirements of library construction and sequencing. BGI does not suggest in using this sample.
- 2. According to BGI's data, the one-time successful rate of library construction is more than 95% for samples of level A and level B.
- 3. According to BGI's data, the risks of library construction for sample of level C or level D are listed below:
- a) The deficiency of the quantity: There may be the risk of failure in library construction and the yield of library of experiment may be too low to sequencing, and the database of low yield for sequencing may lead to poor randomness.
- b) Degradation of sample: It may cause high duplication rate of library and insert fragment will be abnormal."
- c) Pollution by Protein or Insoluble Impurity: It may affect the fragmentation effect, leads to insert size unstable, influence the SNP analysis
- d) RNA contamination: It possibly effects the DNA concentration quantitative accuracy.
- 4. If the partner insists on using the sample of level C or level D, the risk and responsibility is taken by the cooperative partner.
- 5. Other notes:
- a) Some samples are contaminated by proteins.
- b) Some samples are contaminated by RNA.

4. Appendix

Appendix 1: Test results of Qubit Fluorometer or Microplate Reader

Appendix 2: Test results of Agarose Gel Electrophoresis

Appendix 3: Original information of sample

5. Statement

- 1. The results shown in this report refer only to the sample of the report unless otherwise stated.
- 2. This test report cannot be copied partly without the prior written permission of the Lab.

Appendix 1: Test results of Qubit Fluorometer or Microplate Reader

1. Pre-treatment

After the sample melted the ice, centrifuged and fully mixed, take appropriate samples for testing.

2. Test Result

Sample Number	Test Instrument	Test Kit	Dilution Ratio(×)	Test Volume (μL)	$\begin{array}{c} \text{Test} \\ \text{Concentration(ng/} \\ \mu L) \end{array}$	Concentration of original sample(ng/µL)	Remark
8521512001643	Qubit	DNA BR	1	1	154	154	
8521512001643	Qubit	DNA BR	1	1	154	154	
8521512001643	Qubit	DNA BR	1	1	154	154	
8521512001643	Qubit	DNA BR	1	1	154	154	
8521512001643	Qubit	DNA BR	1	1	154	154	
8521512001644	Qubit	DNA BR	1	1	428	428	
8521512001644	Qubit	DNA BR	1	1	428	428	
	Number 8521512001643 8521512001643 8521512001643 8521512001643 8521512001644	Number Instrument 8521512001643 Qubit 8521512001643 Qubit 8521512001643 Qubit 8521512001643 Qubit 8521512001643 Qubit 8521512001644 Qubit	Number Instrument Kit 8521512001643 Qubit DNA BR 8521512001644 Qubit DNA BR 8521512001644 Qubit DNA BR 8521512001644 Qubit DNA BR 8521512001644 Qubit DNA	Number Instrument Kit Ratio(×) 8521512001643 Qubit DNA BR 1 8521512001644 Qubit DNA BR 1 8521512001644 Qubit DNA BR 1 8521512001644 Qubit DNA BR 1	Number Instrument Kit Ratio(×) (μ L) 8521512001643 Qubit DNA BR 1 1 8521512001644 Qubit DNA BR 1 1 8521512001644 Qubit DNA BR 1 1 8521512001644 Qubit DNA BR 1 1	Sample Test Instrument Kit Ratio(×) Test Olime Concentration(ng/μL)	Sample Test Test Dilution Ratio (×) Concentration (ng/ μL) Concentration of original sample (ng/ μL)

oyster_02	8521512001644	Qubit	DNA BR	1	1	428	428	
oyster_02	8521512001644	Qubit	DNA BR	1	1	428	428	
oyster_02	8521512001644	Qubit	DNA BR	1	1	428	428	
oyster_03	8521512001645	Qubit	DNA BR	1	1	296	296	
oyster_03	8521512001645	Qubit	DNA BR	1	1	296	296	
oyster_03	8521512001645	Qubit	DNA BR	1	1	296	296	
oyster_03	8521512001645	Qubit	DNA BR	1	1	296	296	
oyster_03	8521512001645	Qubit	DNA BR	1	1	296	296	

Appendix 2: Test results of agarose gel electrophoresis

1. Pre-treatment

After the sample melted the ice, centrifuged and fully mixed, take appropriate samples for testing.


2. Test Parameter


Concentration of Agarose Gel: 1 %; Voltage: 150 V; Electrophoresis

Time: 40 min

3. Test Result

(1) Electrophoretogram:





Lane No.	Sample Name	Dilution Ratio(\times)	Test Volume(µL)	Sample Integrity	Remark
M1	λ-Hind III digest(Takara)	1	3		
1	oyster_01	1	0.65	Degraded Moderate	RNA contamination
1	oyster_01	1	0.65	Degraded slightly	
1	oyster_01	1	0.65	Degraded slightly	
1	oyster_01	1	0.65	Degraded slightly	
1	oyster_01	1	0.65	Degraded slightly	
2	oyster_02	10	2. 34	Degraded slightly	RNA contamination
2	oyster_02	10	2. 34	Degraded slightly	
2	oyster_02	10	2. 34	Degraded slightly	
2	oyster_02	10	2. 34	Degraded Moderate	
2	oyster_02	10	2. 34	Degraded slightly	
3	oyster_03	1	0. 34	Degraded slightly	RNA contamination
3	oyster_03	1	0. 34	Degraded slightly	RNA contamination
3	oyster_03	1	0. 34	Degraded slightly	RNA contamination
3	oyster_03	1	0.34	Degraded Moderate	RNA contamination
3	oyster_03	1	0. 34	Degraded slightly	RNA contamination
M2	D2000 (Tiangen)	1	6		

Appendix 3: Original information of sample

Sample Typ	Sample Type:										
Genome DNA											
Sample status:											
Dissolved in 10mM Tris-HCl											
Further Information:											
Sample Name	Species	No. of Tubes	Concentration(ng/ μL)	Volume(μL)	Total Quantity(μg)	Fragment Size	OD260/280	OD260/230	Remark		
oyster_01	Ostrea lurida	1	137	200	27. 4		0	0			
oyster_01	Ostrea lurida	1	137	200	27. 4		0	0			
oyster_01	Ostrea lurida	1	137	200	27. 4		0	0			
oyster_01	Ostrea lurida	1	137	200	27. 4		0	0			
oyster_01	Ostrea lurida	1	137	200	27. 4		0	0			
oyster_02	Ostrea lurida	1	74. 3	100	7. 4		0	0			
oyster_02	Ostrea lurida	1	74. 3	100	7. 4		0	0			
oyster_02	Ostrea lurida	1	74. 3	100	7. 4		0	0			
oyster_02	Ostrea lurida	1	74. 3	100	7. 4		0	0			
oyster_02	Ostrea lurida	1	74. 3	100	7. 4		0	0			
oyster_03	Ostrea lurida	1	219	75	16. 4		0	0			
oyster_03	Ostrea lurida	1	219	75	16. 4		0	0			
oyster_03	Ostrea lurida	1	219	75	16. 4		0	0			
oyster_03	Ostrea lurida	1	219	75	16. 4		0	0			
oyster_03	Ostrea lurida	1	219	75	16. 4		0	0			

Report End